



Red Cell Dosage Calculator

User Guide



S. Grey, Version 2 (Dec 2020) Ref: BT-SW-17

... for a **better** Bolton

Introduction

- The empiric nature of red cell dosing in non-bleeding patients with normovolaemic anaemia can lead to under or over-transfusion when aiming to meet a post-transfusion haemoglobin target, especially when the patient's body weight is not accounted for.
- It is always important to clinically re-evaluate patients after a single unit of red cells regardless of the total number administered. However repeat haemoglobin testing between units has a time and resource impact. This could be avoided if the volume required to meet the target haemoglobin level could be reliably predicted in the context of body weight.
- A web app was developed, technically validated and **CE marked as a class 1 medical device** in 2016 by Bolton NHS Foundation Trust, and clinically implemented in 2017 (version 1), and redeveloped in 2019 which would allow sharing with other NHS organisations (version 2).
- The implementation strategy was supported at organisational executive level and included hosting the App on the Trust intranet, promoting individual desktop access, clinical and laboratory education and awareness, embedding in blood transfusion policy, redesign of the transfusion red cell requesting process, and monthly quality impact assessment audit with staff support and feedback.

Computer says no?


- The calculator may well sometimes be contrary to the clinician's view but the calculator is intended to give the best guidance to prescribe sufficient red cells to meet a target Hb. We know that patients may be symptomatic on the margins of the transfusion triggers (70g/L or 80g/L) which is why a 5g/L 'tolerance' is built into the calculator to avoid this.
- Importantly, the Red Cell Dosage Calculator is **Decision-Support Software**. There may be a good clinical reason to prescribe a different volume: for example if two units is the calculated dose but you feel the patient is at risk of circulatory overload and only want one unit, this is clinically justified and can be documented.
- A **Transfusion-Associated Circulatory Overload (TACO)** risk assessment should be performed before every transfusion episode, and this should also be part of the decision on the appropriate red cell volume to prescribe.

Is it effective?

- After the first month of implementation 56% of red cell issues had evidence of the App being used to calculate the dose, rising to 88% by month nine. There were 581 red cell transfusions where patients received the same number of units as calculated by the web app. One hundred and forty five were excluded from the data because there was either no post-transfusion haemoglobin value, were found to be bleeding on retrospective audit, had profound anaemia and would have required more than two units to reach their target haemoglobin, or received fewer units than calculated and issued.
- Of the remaining 436 cases, 388 (89%) reached their post-transfusion target haemoglobin level. Twenty-four cases (5.5%) did not meet the target, and 24 (5.5%) exceeded the target. At financial month 8 (2017/18), 3016 units had been transfused in total with a prediction that 3348 would be transfused for normovolaemic anaemia by financial year end, representing a projected 16.5% reduction in red cell usage (18% actual).


Intranet Access

The App can be hosted on an organisation's intranet page



... for a **better** Bolton

Personal Desktop Access



Enter the home page of the calculator and copy the url address

Right click on your desktop
Click 'New'...'Shortcut'
Paste the url address

Click 'next'

Name the shortcut:
'Red Cell Dosage Calculator'

Click 'finish'

Exclude Non-Appropriate Patients

To be used on all non-bleeding adult patients who are not transfusion-dependent

The screenshot shows the 'Red Cell Dosage Calculator' web application. At the top left is the logo and title 'Red Cell Dosage Calculator'. At the top right is a 'Contact Us' button. Below the title bar is a breadcrumb trail 'Home / Pre-Step' and the text 'Software Version: 2.0'. The main content area contains a question: 'Is The Patient:' followed by a bulleted list: 'Actively bleeding?', 'On a regular transfusion programme?', and 'Paediatric/neonatal?'. Below the list are two buttons: a red 'Yes' button and a green 'No' button. At the bottom of the page, there are links for 'Privacy Policy', 'Cookie Policy', 'Accessibility', and 'Terms and Conditions', along with the copyright notice '©2020 Technical Solutions Worldwide Ltd.'. On the right side of the footer, there is a CE mark, the NHS Bolton logo, and the text 'Class I Medical Device (Decision Support Software)'.

Enter Details...

The screenshot shows the 'Red Cell Dosage Calculator' interface. At the top left is the calculator's logo, and at the top right is a 'Contact Us' button. Below the header is a breadcrumb trail: Home / Pre-Step / Calculator / Result. The main content area is divided into sections: 'Pre-Transfusion Haemoglobin and Body Weight' with input fields for haemoglobin (68 g/L) and body weight (61 Kg); 'Cause of Anaemia' with a text field containing 'Anaemia caused by surgical blood loss, or GI / Gynae / Obstetric bleeding (without current active bleeding)'; and 'Co-morbidities and Risk Factors' with a text field containing 'Risk factors for chronic cardio-vascular disease (e.g. > 70 yrs old, peripheral vascular disease, diabetes)'. A 'Results' section displays a green box with the calculated dose (293 mLs) and the required number of units (1). A 'New Calculation' button is located at the bottom of the main content area. The footer contains links for Privacy Policy, Cookie Policy, Accessibility, and Terms and Conditions, along with copyright information for Technical Solutions Worldwide Ltd. Regulatory logos for CE and NHS Bolton are also present, along with the text 'Class I Medical Device (Decision Support Software)'.

... for a **better** Bolton

Example: Borderline Pre- Transfusion Hb

The screenshot displays the 'Red Cell Dosage Calculator' interface. At the top left is the logo and title 'Red Cell Dosage Calculator', and at the top right is a 'Contact Us' button. Below the title is the version number 'Version: 2.0'. A breadcrumb trail reads 'Home / Pre-Step / Calculator / Result'. The main content area is divided into sections: 'Pre-Transfusion Haemoglobin and Body Weight' with input fields for 'The patient's pre-transfusion haemoglobin: 72 g/L' and 'The patient's body weight: 69 Kg'; 'Cause of Anaemia' with the text 'Anaemia caused by surgical blood loss, or GI / Gynae / Obstetric bleeding (without current active bleeding)'; and 'Co-morbidities and Risk Factors' with the text 'There are no co-morbidities or risk factors for this patient.' Below this is a 'Results' section containing a green box with the text: 'Transfusion not indicated as the current haemoglobin equals or exceeds the target haemoglobin for the indication given. However, if the patient is significantly symptomatic, consider one unit.' At the bottom of the main content area is a 'New Calculation' button. The footer contains links for 'Privacy Policy', 'Cookie Policy', 'Accessibility', and 'Terms and Conditions', along with the copyright notice '©2021 Technical Solutions Worldwide Ltd.', the CE mark, the NHS Bolton logo, and the text 'NHS Foundation Trust' and 'Class I Medical Device (Decision Support Software)'.

... for a better Bolton

**Example: >1
unit still
requires
review after
1st unit**

Red Cell Dosage Calculator [Contact Us](#)

Version: 2.0

[Home](#) / [Pre-Step](#) / [Calculator](#) / [Result](#)

Pre-Transfusion Haemoglobin and Body Weight

The patient's pre-transfusion haemoglobin: 62 g/L

The patient's body weight: 95 Kg

Cause of Anaemia

Anaemia caused by surgical blood loss, or GI / Gynae / Obstetric bleeding (without current active bleeding).

Co-morbidities and Risk Factors

Risk factors for chronic cardio-vascular disease (e.g. > 70 yrs old, peripheral vascular disease, diabetes).

Results

Calculated Dose = 684 mLs.

Required Number of Units of Red Cells = 2

Transfuse 1 unit then re-evaluate the symptoms / haemoglobin before administering the second unit. This is under the calculated dose.

Medical management of anaemia is recommended if the target Hb is not met following transfusion.

[New Calculation](#)

[Privacy Policy](#) [Cookie Policy](#) [Accessibility](#) [Terms and Conditions](#)

©2021 Technical Solutions Worldwide Ltd.

Class I Medical Device (Decision Support Software)

... for a better Bolton

Example: Low body weight patients

Red Cell Dosage Calculator Contact Us

Version: 2.0

[Home](#) / [Pre-Step](#) / [Calculator](#) / [Result](#)

Pre-Transfusion Haemoglobin and Body Weight

The patient's pre-transfusion haemoglobin: 73 g/L

The patient's body weight: 44 Kg

Cause of Anaemia

Chronic anaemia that is not medically reversible.
Anaemia caused by chemotherapy.
Anaemia caused by bone marrow failure.

Co-morbidities and Risk Factors

Risk factors for chronic cardio-vascular disease (e.g. > 70 yrs old, peripheral vascular disease, diabetes).

Results

Calculated Dose = 123 mLs.

This is less than the volume of a single unit of red cells. Consider whether transfusion is necessary.
If the patient is < 50Kg, consider prescribing in mLs rather than units (and administer via a pump).

[New Calculation](#)

[Privacy Policy](#) [Cookie Policy](#) [Accessibility](#) [Terms and Conditions](#)
©2021 Technical Solutions Worldwide Ltd.

Class I Medical Device (Decision Support Software)

Comments and Help

Clinical Comments/Questions

(non-technical or contractual)

htt@boltonft.nhs.uk

Technical and Contractual support:

Contact Technical Solutions Worldwide Ltd via the app
website using the web form

www.rcdcalculator.co.uk

Project Lead:

Dr Sharran Grey (Honorary Haematology
Consultant Clinical Scientist)

Sharran.grey@boltonft.nhs.uk